

Speed Dial / ServiceLINK

When pressed, the button initiates a call to a pre-programmed phone number like, loved ones, friends or doctor's office.

Check Button

Press to hear time of day. Also turns off and on reminders

Ultra Sensitive Microphone

Lets us hear you from anywhere in your home



Pendants

Our **SMART PENDANTS** can be pressed to answer an incoming phone call or disconnect at the end of a phone conversation.



Speaker

This speaker is what makes it possible for you to hear our operators in other rooms of the home

Call Button

Pressing this button, puts you in contact with the monitoring center on a non-emergency basis